



6 August 2004

Complete solution for stability testing management

Autoscribe announces the launch of Matrix Stability Express TK, a new software package designed in conjunction with several pharmaceutical companies to manage and report on protocol driven studies associated with stability or shelf-life testing.

The protocol designer allows the entry and subsequent management of all study data. Furthermore, information is stored in one location for easy access by authorised users. Stability Express also manages storage room operations for quick and easy retrieval of samples and the cycling of conditions, where required. An inventory management capability ensures sufficient material is available for an entire study and prevents undue wastage through the excess storage of containers. As samples are taken from a store room, information is automatically registered onto the results system and the appropriate tests and limits allocated according to the protocol requirements. This saves significant operator time and raises the efficiency of stability testing exercises.

Matrix Stability Express TK comes with a set of standard reports, together with a link to JMP statistical analysis software for calculating shelf life projections, accelerated shelf life and more. Audit trail generation and other features assist in the compliance with FDA regulatory requirements such as 21 CFR Part 11.

Autoscribe is offering a choice of two validation service plans. The first is a set of tools that helps users validate the system using prepared documentation. The second is a full validation service that involves visits to the user site for the execution of validation scripts.

Matrix Stability Express TK provides everything from protocol design and approval to completed trend plots and reports.

About Autoscribe

Autoscribe Limited is a world leader in the development and supply of successful software solutions for LIMS (Laboratory Information Management Systems) and the scientific laboratory and business markets.

Autoscribe pioneered the use of innovative configuration tools within the LIMS market and today our Matrix Plus LIMS software is seen as the industry standard to which others aspire.

Autoscribe offers a wide portfolio of products and services including systems for LIMS, sample tracking, instrument calibration and maintenance, stability testing, customer complaints management, ISO9000 protocol driven training records management and resource planning. Our customers include leading analytical, oil, pharmaceutical, biotechnology, chemical, and consumer products companies as well as academic research institutes, hospitals and many others.

Founded in 1981, Autoscribe is a privately owned UK company. Worldwide we are represented by a network of business partners in regions including the USA, France, Benelux, Germany, Switzerland, Scandinavia and Australia. Our organisation includes a dedicated team of developers, technical support specialists and sales & marketing professionals.

Further Information

For further information please contact us at:

Autoscribe UK Ltd.
Wellington House
Riseley Business Park
Basingstoke Road
Riseley
Reading
Berkshire
RG7 1NW
UK

Telephone +44 (0) 118 984 0610
Fax +44 (0) 118 984 0611

Email: matrix@autoscribe.co.uk

Or visit our website at <http://www.autoscribe.co.uk> and complete the information request form