

| <b>LIMS Terminology</b>  |  |
|--|--|
| <b>LIMS</b>  | Laboratory Information Management System   |
| <b>CDS</b>   | Chromatography Data System   |
| <b>PREPARATION</b>   | A step after receiving a sample into the lab when the sample is processed (maybe extracted) to be ready for analysis.  |
| <b>SAMPLE LIFECYCLE</b>  | The lifecycle / workflow whereby samples, batches, lots are registered, prepared, received, tested, validated, approved and final reports / certificates of analysis are produced. |
| <b>SUBSTANCE</b>   | Product, Material  |
| <b>REGISTRATION</b>  | Log, Logging, Log Sample   |
| <b>Regulatory Terminology</b>                                  |  |
| <b>21 CFR Part 11,<br/>Part 11, E-Records,<br/>E-Signature</b> | Federal Food and Drug Administration regulation 21CFR Part 11, governing electronic records and electronic signatures.   |
| <b>GMP</b>   | Good Manufacturing Practice  |
| <b>SOP / OP</b>  | Standard Operating Procedure, Operating Procedure  |
| <b>IQ</b>  | Installation Qualification   |
| <b>OQ</b>  | Operational Qualification  |
| <b>URS</b>   | User Requirements Specification  |
| <b>FRS</b>   | Functional Requirements Specification  |
| <b>SDS</b>   | System Design Specification  |
| <b>SATC</b>  | Site acceptance test checklist   |
| <b>Computer / Networking Terminology</b>                       |  |
| <b>IP</b>  | Internet Protocol  |
| <b>TCP/IP</b>  | Transmission Control Protocol/Internet Protocol, the communications protocols used to connect hosts on the Internet. The de facto standard for transmitting data over networks.    |
| <b>OS</b>  | Operating System   |