

Autoscribe News Release – PR0504

FOR IMMEDIATE RELEASE
5th April 2005, Riseley, UK

Manage your resources more effectively with a new version of Resource Planner from Autoscribe UK Ltd

Effective deployment of valuable resources is fundamental to the success of any business. To help businesses raise productivity through better resource tracking and planning, Autoscribe is pleased to announce that a new version of its Resource Planner software is now available.

Built using the latest Microsoft.Net framework, Resource Planner has the look and feel of Matrix LIMS, Autoscribe's primary laboratory management system. A unique aspect of the system is the use of artificial intelligence algorithms to generate plans for maximum performance and results.

Resource Planner looks at any new or outstanding tasks and suggests an optimal plan for completion based on workload priority, task completion times, and resource availability. Short term planning is achieved through automatic weekly scheduling and work sheet generation, while longer term planning is possible through due date checking and task completion time allocation.

Once relevant resource information has been entered into Resource Planner, a plan can be generated using the scheduling algorithm. The generated plan can then be activated so the resources in the plan can be booked and subsequent project plans generated. The output from Resource Planner includes the generated plans and reports, as well as an expertise tracker and equipment/instrument tracker. Information is presented in clear, highly intuitive graphs for quick and easy assessment.

Often resources available can change at short notice, which can result in delays in the completion of important and time-sensitive projects. To address this issue Resource Planner uses 'what-if' scenarios by recalculating work completion times, and adjusting the priority of tasks accordingly.

Resource Planner is available as a stand alone product or can be integrated as a module of Matrix LIMS, an existing LIMS or other Information Management System.

EDITORIAL CONTACTS:

John Boother, Managing Director
E-mail: matrix@autoscribe.co.uk

Autoscribe - The Winning Formula!

About Autoscribe

Autoscribe Limited is a world leader in the development and supply of successful software solutions for LIMS (Laboratory Information Management Systems) and the scientific laboratory and business markets.

Autoscribe pioneered the use of innovative configuration tools within the LIMS market and today our Matrix Plus LIMS software is seen as the industry standard to which others aspire.

Autoscribe offers a wide portfolio of products and services including systems for LIMS, sample tracking, instrument calibration and maintenance, stability testing, customer complaints management, ISO9000 protocol driven training records management and resource planning. Our customers include leading analytical, oil, pharmaceutical, biotechnology, chemical, and consumer products companies as well as academic research institutes, hospitals and many others.

Founded in 1981, Autoscribe is a privately owned UK company. Worldwide we are represented by a network of business partners in regions including the USA, France, Benelux, Germany, Switzerland, Scandinavia and Australia. Our organisation includes a dedicated team of developers, technical support specialists and sales & marketing professionals.

For more information, see Autoscribe on the worldwide web at www.autoscribe.co.uk

Further Enquiries

Autoscribe UK Ltd.

Wellington House
Riseley Business Park
Basingstoke Road
Riseley
Reading
Berkshire
RG7 1NW
UK

Telephone: +44 (0) 118 984 0610

Fax: +44 (0) 118 984 0611

Email: matrix@autoscribe.co.uk

Alternatively, please complete our information request form at:
<http://www.autoscribe.co.uk/requestforinfo.shtml>

Reader Enquiries

Please send reader enquiries from this news release directly to matrix@autoscribe.co.uk